

```

*****
*DSYNSM5 - DESY MEDPAR OUTPUT COPYBOOK WITH 50 BYTE HEADER
*          LRECL=1424
*          THIS FORMAT IS FOR BOTH IDENTIFIABLE AND LDS FILES.
*****
*****01  MEDPAR-1424-NCH-REC.
          02  MEDPAR-HEADER-GRP.
              03  MEDPAR-HEADER-VALUE                PIC X(30) .
              03  MEDPAR-HEADER-TBD                  PIC X(11) .
              03  MEDPAR-HEADER-LINK                  PIC X(09) .
          02  MEDPAR-NCH-CLM-TYPE-CD                  PIC X(02) .
          02  MEDPAR-CLM-LCTR-NUM-GRP.
              03  MEDPAR-BENE-CLM-ACNT-NUM            PIC X(09) .
              03  MEDPAR-CTGRY-EQTBL-BIC-CD          PIC X(02) .
          02  MEDPAR-BENE-AGE-CNT                      PIC 9(03) .
          02  MEDPAR-BENE-SEX-CD                      PIC X(01) .
          02  MEDPAR-BENE-RACE-CD                    PIC X(01) .
          02  MEDPAR-BENE-MDCR-STUS-CD                PIC X(02) .
          02  MEDPAR-BENE-RSDNC-SSA-STATE-CD          PIC X(02) .
          02  MEDPAR-BENE-RSDNC-SSA-CNTY-CD          PIC X(03) .
          02  MEDPAR-BENE-MLG-CNTCT-ZIP-CD           PIC X(05) .
          02  MEDPAR-ADMSN-DAY-CD                    PIC 9(01) .
          02  MEDPAR-BENE-DSCHRG-STUS-CD              PIC X(01) .
          02  MEDPAR-GHO-PD-CD                      PIC X(01) .
          02  MEDPAR-PPS-IND-CD                      PIC X(01) .
NEW      02  MEDPAR-ORG-NPI-NUM                      PIC X(10) .
          02  MEDPAR-PRVDR-NUM-GRP.
              03  MEDPAR-PRVDR-STATE-CD              PIC 9(02) .
              03  MEDPAR-PRVDR-NUM-3RD-CD            PIC X(01) .
              03  MEDPAR-PRVDR-NUM-SRL-CD            PIC X(03) .
          02  MEDPAR-PRVDR-NUM-SPCL-UNIT-CD          PIC X(01) .
          02  MEDPAR-SS-LS-SNF-IND-CD                PIC X(01) .
          02  MEDPAR-STAY-FINL-ACTN-CLM-CNT          PIC +9(03) .
          02  MEDPAR-LTST-CLM-ACRTN-DT              PIC +9(07) .
          02  MEDPAR-BENE-MDCR-BNFT-EXHST-DT         PIC +9(07) .
          02  MEDPAR-SNF-QUALN-FROM-DT               PIC +9(07) .
          02  MEDPAR-SNF-QUALN-THRU-DT               PIC +9(07) .
          02  MEDPAR-ADMSN-DT                       PIC +9(07) .
          02  MEDPAR-DSCHRG-DT                      PIC +9(07) .
          02  MEDPAR-CVR-LVL-CARE-THRU-DT            PIC +9(07) .
          02  MEDPAR-BENE-DEATH-DT                  PIC +9(07) .
          02  MEDPAR-BENE-DEATH-DT-VRFY-CD          PIC X(01) .
          02  MEDPAR-INTRNL-USE-SSI-GRP.
              03  MEDPAR-INTRNL-USE-SSI-IND-CD        PIC X(01) .
UPDATE   03  MEDPAR-INTRNL-USE-SSI-DAY-CNT          PIC +9(05) .
              03  MEDPAR-INTRNL-USE-SSI-DATA          PIC X(01) .
          02  MEDPAR-LOS-DAY-CNT                    PIC +9(05) .
          02  MEDPAR-OUTLIER-DAY-CNT                 PIC +9(03) .
UPDATE   02  MEDPAR-UTLZTN-DAY-CNT                 PIC +9(05) .
          02  MEDPAR-TOT-COINSRNC-DAY-CNT           PIC +9(03) .
          02  MEDPAR-BENE-LRD-USE-CNT               PIC +9(03) .
          02  MEDPAR-FILLER-3                      PIC X(12) .
          02  MEDPAR-BENE-PTA-COINSRNC-AMT          PIC +9(07) .
          02  MEDPAR-BENE-IP-DDCTBL-AMT            PIC +9(07) .

```

02	MEDPAR-BENE-BLOOD-DDCTBL-AMT	PIC +9(07).
02	MEDPAR-BENE-PRMRY-PYR-AMT	PIC +9(07).
02	MEDPAR-DRG-OUTLIER-PMT-AMT	PIC +9(07).
02	MEDPAR-IP-DSRPRTNT-SHR-AMT	PIC +9(07).
02	MEDPAR-IME-AMT	PIC +9(07).
02	MEDPAR-DRG-PRICE-AMT	PIC +9(07).
02	MEDPAR-PASS-THRU-AMT	PIC +9(07).
02	MEDPAR-TOT-PPS-CPTL-AMT	PIC +9(07).
02	MEDPAR-IP-LOW-VOL-PMT-AMT	PIC +9(7)V99.
02	MEDPAR-FILLER-4	PIC X(14).
02	MEDPAR-TOT-CHRG-AMT	PIC +9(07).
02	MEDPAR-TOT-CVR-CHRG-AMT	PIC +9(07).
02	MEDPAR-MDCR-PMT-AMT	PIC +9(07).
02	MEDPAR-ACMDTNS-TOT-CHRG-AMT	PIC +9(07).
02	MEDPAR-DPRTMNTL-TOT-CHRG-AMT	PIC +9(07).
02	MEDPAR-ACMDTNS-DAYS-GRP.	
03	MEDPAR-PRVT-ROOM-DAY-CNT	PIC +9(03).
03	MEDPAR-SEMI PRVT-ROOM-DAY-CNT	PIC +9(03).
03	MEDPAR-WARD-DAY-CNT	PIC +9(03).
03	MEDPAR-INTNSV-CARE-DAY-CNT	PIC +9(03).
03	MEDPAR-CRNRV-CARE-DAY-CNT	PIC +9(03).
02	MEDPAR-ACMDTNS-CHRG-GRP.	
03	MEDPAR-PRVT-ROOM-CHRG-AMT	PIC +9(07).
03	MEDPAR-SEMI PRVT-ROOM-CHRG-AMT	PIC +9(07).
03	MEDPAR-WARD-CHRG-AMT	PIC +9(07).
03	MEDPAR-INTNSV-CARE-CHRG-AMT	PIC +9(07).
03	MEDPAR-CRNRV-CARE-CHRG-AMT	PIC +9(07).
02	MEDPAR-SRVC-CHRG-GRP.	
03	MEDPAR-OTHR-SRVC-CHRG-AMT	PIC +9(07).
03	MEDPAR-PHRMCY-CHRG-AMT	PIC +9(07).
03	MEDPAR-MDCL-SUPPLY-CHRG-AMT	PIC +9(07).
03	MEDPAR-DME-CHRG-AMT	PIC +9(07).
03	MEDPAR-USED-DME-CHRG-AMT	PIC +9(07).
03	MEDPAR-PHYS-THRPY-CHRG-AMT	PIC +9(07).
03	MEDPAR-OCPTNL-THRPY-CHRG-AMT	PIC +9(07).
03	MEDPAR-SPCH-PTHLGY-CHRG-AMT	PIC +9(07).
03	MEDPAR-INHLTN-THRPY-CHRG-AMT	PIC +9(07).
03	MEDPAR-BLOOD-CHRG-AMT	PIC +9(07).
03	MEDPAR-BLOOD-ADMIN-CHRG-AMT	PIC +9(07).
03	MEDPAR-OPRTG-ROOM-CHRG-AMT	PIC +9(07).
03	MEDPAR-LTHTRPSY-CHRG-AMT	PIC +9(07).
03	MEDPAR-CRDLGY-CHRG-AMT	PIC +9(07).
03	MEDPAR-ANSTHSA-CHRG-AMT	PIC +9(07).
03	MEDPAR-LAB-CHRG-AMT	PIC +9(07).
03	MEDPAR-RDLGY-CHRG-AMT	PIC +9(07).
03	MEDPAR-MRI-CHRG-AMT	PIC +9(07).
03	MEDPAR-OP-SRVC-CHRG-AMT	PIC +9(07).
03	MEDPAR-ER-CHRG-AMT	PIC +9(07).
03	MEDPAR-AMBLNC-CHRG-AMT	PIC +9(07).
03	MEDPAR-PROFNL-FEES-CHRG-AMT	PIC +9(07).
03	MEDPAR-ORGN-ACQSTN-CHRG-AMT	PIC +9(07).
03	MEDPAR-ESRD-REV-SETG-CHRG-AMT	PIC +9(07).
03	MEDPAR-CLNC-VISIT-CHRG-AMT	PIC +9(07).
02	MEDPAR-ACMDTNS-SRVC-IND-GRP.	

03	MEDPAR-ICU-IND-CD	PIC X(01).
03	MEDPAR-CRNR-CARE-IND-CD	PIC X(01).
03	MEDPAR-PHRMCY-IND-CD	PIC 9(01).
03	MEDPAR-TRNSPLNT-IND-CD	PIC 9(01).
03	MEDPAR-RDLGY-IND-GRP.	
04	MEDPAR-RDLGY-ONCLGY-IND-SW	PIC 9(01).
04	MEDPAR-RDLGY-DGNSTC-IND-SW	PIC 9(01).
04	MEDPAR-RDLGY-THRPTC-IND-SW	PIC 9(01).
04	MEDPAR-RDLGY-NUCLR-MDCN-IND-SW	PIC 9(01).
04	MEDPAR-RDLGY-CT-SCAN-IND-SW	PIC 9(01).
04	MEDPAR-RDLGY-OTHR-IMGNG-IND-SW	PIC 9(01).
03	MEDPAR-OP-SRVC-IND-CD	PIC 9(01).
03	MEDPAR-ORGN-ACQSTN-IND-CD	PIC X(02).
03	MEDPAR-ESRD-SETG-IND-CD	OCCURS 00005 TIMES
		PIC X(02).
02	MEDPAR-POA-DGNS-CD-GRP.	
03	MEDPAR-POA-DGNS-CD-CNT	PIC 9(02).
03	MEDPAR-POA-DGNS-IND-CD	OCCURS 00025 TIMES
		PIC X(01).
02	MEDPAR-POA-DGNS-E-CD-GRP.	
03	MEDPAR-POA-DGNS-E-CD-CNT	PIC 9(02).
03	MEDPAR-POA-DGNS-E-IND-CD	OCCURS 00012 TIMES
		PIC X(01).
02	MEDPAR-DGNS-CD-GRP.	
03	MEDPAR-DGNS-CD-CNT	PIC 9(02).
03	MEDPAR-DGNS-VRSN-CD	PIC X(01).
03	MEDPAR-DGNS-CD	OCCURS 00025 TIMES
		PIC X(07).
02	MEDPAR-DGNS-E-CD-GRP.	
03	MEDPAR-DGNS-E-CD-CNT	PIC 9(02).
03	MEDPAR-DGNS-E-VRSN-CD	PIC X(01).
03	MEDPAR-DGNS-E-CD	OCCURS 0012 TIMES
		PIC X(07).
02	MEDPAR-SRGCL-PRCDR-IND-SW	PIC X(01).
02	MEDPAR-SRGCL-PRCDR-GRP.	
03	MEDPAR-SRGCL-PRCDR-CD-CNT	PIC 9(02).
03	MEDPAR-SRGCL-PRCDR-DT-CNT	PIC 9(02).
03	MEDPAR-SRGCL-PRCDR-VRSN-CD	PIC X(01).
03	MEDPAR-SRGCL-PRCDR-CD	OCCURS 00025 TIMES
		PIC X(07).
03	MEDPAR-SRGCL-PRCDR-PRFRM-DT	OCCURS 00025 TIMES
		PIC +9(07).
02	MEDPAR-BLOOD-PT-FRNSH-QTY	PIC +9(03).
02	MEDPAR-BENE-IDENT-CD	PIC X(02).
02	MEDPAR-DRG-CD	PIC 9(03).
02	MEDPAR-DSCHRG-DSTNTN-CD	PIC 9(02).
02	MEDPAR-DRG-OUTLIER-STAY-CD	PIC 9(01).
02	MEDPAR-BENE-PRMRY-PYR-CD	PIC X(01).
02	MEDPAR-ESRD-COND-CD	PIC 9(02).
02	MEDPAR-SRC-IP-ADMSN-CD	PIC X(01).
02	MEDPAR-IP-ADMSN-TYPE-CD	PIC X(01).
02	MEDPAR-FICARR-IDENT-NUM	PIC X(05).
02	MEDPAR-ADMTG-DGNS-CD-GRP.	
03	MEDPAR-ADMTG-DGNS-VRSN-CD	PIC X(01).

03	MEDPAR-ADMTG-DGNS-CD	PIC X(07).
02	MEDPAR-ADMSN-DEATH-DAY-CNT	PIC +9(05).
02	MEDPAR-FILLER-6	PIC X(04).
02	MEDPAR-INTRNL-USE-IPSB-CD	PIC 9(03).
02	MEDPAR-INTRNL-USE-FIL-DT-CD	PIC 9(01).
02	MEDPAR-INTRNL-USE-SMPL-SIZE-CD	PIC 9(01).
02	MEDPAR-WRNG-IND-CD	PIC +9(17).
02	FILLER	PIC X(19).
02	MEDPAR-ORIG-CLM-LCTR-NUM-GRP	PIC X(11).
02	MEDPAR-ACTIVE-XREF-IND	PIC X.
02	MEDPAR-SELECTION-REASON	PIC X.
*****	END OF DSYNSM5	*****

—